



FASEROPTIK


RFQ-Medizintechnik GmbH & Co. KG

Bruderhofstr. 10 - 12
D-78532 Tuttlingen

Telefax (07461) 9617-20
Telefon (07461) 9617-0

E-Mail: info@rfq.de
<http://www.kaltlichtkabel.de>

COMPARISON LIST FIBER LIGHT GUIDES



LINVATEC CONMED

SYSTEM STANDARD

Fibre Bundle Diameter	Length mm	RfQ-Article Cable	Projector Fitting	Endoscope Fitting	Linvatec-Article
3,5	3000	00-21110-63	MORE POSSIBILITIES (SEE BELOW)		C3279
4,9	3000	00-21210-63		C3278	

SYSTEM FUSED

TO BE USED WITH HIGH TEMPERATURE LIGHT SOURCES

Fibre Bundle Diameter	Length mm	RfQ-Article Cable	Projector Fitting	Endoscope Fitting	Linvatec-Article
3,5	3000	00-22330-50	MORE POSSIBILITIES (SEE BELOW)		C3279
4,9	3000	00-22230-50		C3278	

DIFFERENT FITTINGS

LE ACMI	00-21110-66	7456
LE WOLF, DYONICS	00-21110-74	7455
LE STORZ	00-21110-70	7457
LE OLYMPUS	00-21116-62	7458
LA OLYMPUS	00-21010-71	7453
LA WOLF, DYONICS	00-21010-76	7451
LA ACMI	00-21010-67	7452

Fiber Light Guides are not 100% like Cables from Original Manufacturer

Errors excepted

LE=LIGHT INPUT
LA=LIGHT OUTPUT